# Dentimplant

## CERRAHİ KULLANIM KILAVUZU

### DENTİMLANT (DENTCON) SİSTEMİ

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## SURGICAL INSTRUCTION MANUAL

### DENTIMPLANT (DENTCON) SYSTEM

It contains information for the use of fixtures, prosthetic components and surgical instruments.

### GENERAL DESCRIPTION

Implant is a dental screw made of titanium metal intended to be surgically inserted into the bone of the upper or lower jaw arches. The fixture is made of pure titanium metal and supplied sterile.

Shelflife is 4 years.  
For Single Use Only.

### Using DENTIMPLANT (DENTCON) system:

#### Preparations before use

1. The dentist should be fully aware of the product and procedure and should inform the patient about the implant application before use. The patient should be aware of the functional and aesthetic limitations of the implant.
2. Selection and placement of the appropriate implant is an important factor affecting the life of the implant, therefore its indications, contraindications, precautions and warnings should be strictly followed.
3. Before the procedure, the patient's condition should be thoroughly examined and a suitable method should be determined for the patient. Failure in this regard may result in implant loss.
4. Examine the x-ray film and choose the right product for the patient.
5. Check the product package status and examine the expiry date information and package damage status.
6. The product is presented to you as "sterile", do not use products with opened or damaged packaging.

#### Packaging

Packaging Fixtures are packaged in tubes with titanium holders and sent to the gamma irradiation facility for sterilization. Products returning from sterilization are packaged in final boxes.

#### Labeling

Fixtures come with tube labels, double labels, and box labels. These labels allow for three-step lot number

control and provide implant diameter and length information.

#### Caution

- All models in DENTIMPLANT (DENTCON) system must be applied according to DENTIMPLANT (DENTCON) Surgery and Prosthesis protocols. It is recommended to wait at least 3 months (12 weeks) for osseointegration after the surgical procedure.
- Implants are in sterile packaging.
- Do not open the cap of the tube in any way until the moment of use. Never use tubes that have been opened or deformed.
- Keep out of the reach of children.
- Implants should be placed in a suitable diameter, sufficient number and in an axis compatible with the tooth arrangement.
- Use mills suitable for implant diameter and length.
- Inform the patient before and after the surgical operation.
- The health condition of the patient must be appropriate for the operation.
- Against the risk of contamination, the implant should be placed in the slot opened without touching any place as soon as it is removed from the sterile tube.
- In implant treatments, it should be observed that the patient completes the bone development.

**Limitations for Bone and Level Implants Note:** The implant carrier must be used with a maximum speed of 25 rpm/min when used with a hand tool, surgical motor.

1. Insert the implant into the implant bed with a maximum torque of 45 Ncm. Exceeding the recommended torque value may cause bone necrosis or damage/failure of the implant.
2. Failure to follow the drilling protocol during drilling may cause the implant to become stuck without being fully inserted. In this case, use the implant motor in reverse mode or manually rotate it counterclockwise with the ratchet to remove the implant and put it back into the inner package. Replace the implant by completing the drill protocol. Caution: A torque value of at least 35 Ncm must be reached for immediate loading.

#### Warning

Previously used implants cannot be reused. Use only DENTIMPLANT (DENTCON) System parts and surgical kit. The company does not accept any responsibility in case of any use of a part of another system.

#### Traceability

Each implant package has a Lot number on it. In addition, two Lot number labels are included in each implant package. For retrospective traceability of the product, these Lot number labels must be attached to the patient's file and panoramic x-ray.

#### Patient Population

Dental Implants are designed for use in patients with complete or partial edentulism, who have completed their growth and development, do not have the conditions specified in the contraindications. (Patients older than 18 years of age who have completed jawbone development.)

#### Duration of the Benefit

The duration of benefit of an implant can be defined as the time it maintains its function. An implant failure can be caused by screw loosening, fracture, eruption etc. Success rates of dental implants vary, depending on where in the jaw the implants are placed but, in general, dental implants have a success rate of up to 98%. With proper care, implants can last a lifetime (10-30 years).

#### Clinical Claims

##### Implant Survival:

Sustainability of long-term functional rehabilitation is achieved. Implant survival rates are reported as 83%-100% overall and typically 96%-99% in routine cases.

##### Implant Loss:

Continuity of treatment success is high. Implant loss is reported at 0%-3% in routine cases and 8%-22% in complex/atrophic cases.

**Marginal Bone Stability (MBL):** Peri-implant bone support is preserved. Marginal bone loss is reported as 0.02-1.22 mm with a long-term typical value of ~0.4-0.7 mm.

##### Peri-implantitis / Infection:

Biological safety and stability are maintained. Peri-implant infection rates are generally 0%-2.7%.

##### Soft Tissue Stability / Esthetic Outcome:

Esthetic appearance is preserved. >1 mm recession may be reported in up to 11% of cases; average recession is 0.2-1.3 mm.

**Mechanical / Prosthetic Complications:** Continuity of functional restoration is maintained.

Mechanical/prosthetic complications may occur at 1%-12% and are manageable.

##### Prosthesis / Crown Survival:

Sustainability of chewing function is achieved. Prosthesis/crown survival has been reported as 94.7%-100%.

##### Patient Satisfaction / Quality of Life:

Improvement in quality of life is achieved. Patient satisfaction rates range between 85%-97%.

## MARKETING CLINICAL CLAIMS

### High Implant Survival / High Success

Implant survival ranges between 83%-100% overall and is typically 96%-99% in routine cases.

### Low Implant Loss / Low Failure

Implant loss is reported at 0%-3% in routine cases and 8%-22% in complex/atrophic cases.

### Bone Stability Is Preserved

Marginal bone loss ranges between 0.02-1.22 mm with a long-term typical value of approximately 0.4-0.7 mm.

**Low Infection / Peri-implant Disease Rate** Peri-implant infection rates are generally between 0%-2.7%

### Soft Tissue / Esthetic Stability

>1 mm recession may be reported in up to 11% of cases; average soft tissue change is 0.2-1.3 mm.

### Low Mechanical/Prosthetic Problems, Function Continues

Mechanical/prosthetic complications may occur at 1%-12% and are generally manageable with maintenance.

**High Prosthesis / Crown Success** Prosthesis/crown survival has been reported between 94.7%-100%.

### Patient Satisfaction / Quality of Life Improves

Patient satisfaction and quality of life indicators range between 85%-97%

## Indications

DENTIMPLANTS (DENTCON) are medical device products that are surgically placed in the lower and upper jaw bones to support artificial dental prostheses used to restore the chewing function of the patient.

## Contraindications

High blood pressure, cardiovascular diseases, diabetes, bone metabolism disorders, uncontrolled bleeding diseases, insufficient wound healing capacity, inadequate oral hygiene, serious internal medical problems, incomplete development in the lower and upper jaw, poor general health, non-cooperative and non-motivated patient, drug, alcohol or cigarette addiction, psychosis, long-term treatment-resistant functional disorders, xerostomia (dry mouth), granulocytopenia, EhlersDanlos syndrome, radiation therapy, osteoradionecrosis, weak immune system, kidney failure, organ transplantation, fibrous dysplasia, Crohn's disease, use of steroids, corticosteroids or anticonvulsants, prophylactic antibiotics, creatinine, serum calcium, titanium allergy, uncontrolled endocrine disorders, anticoagulation drugs, hemorrhagic diathesis, bruxism, parafunctional habits, unfavorable anatomical bone conditions, odontitis, temporomandibular joint disorders, treatable pathological diseases of the jaw or oral mucosa changes, pregnancy, breastfeeding, osteoporotic fracture, respiratory disease, thyroid or parathyroid diseases. Patients with active osteolytic, inflammation, diagnosed malignant, nodular growths, tenderness, unexplained swelling of the head or neck, infectious condition at the implant part. Unrealistic patient expectations. Inaccessible prosthetic reconstructions. Inadequate training of the practitioner.

## Risks

Perforation of the nasal and maxillary sinuses, local and systemic infections, soft tissue perforations, nerve damage. Temporary pain and swelling, speech disorders and gingivitis may occur due to implant placement.

Long-term problems include nerve, local or systemic bacterial infections, and infectious diseases in susceptible individuals. The existing tooth alignment may be compromised as a result of improper implant placement.

## Frequency of Occurrence of the Most Common Identified Risks

**Implant Loss / Osseointegration Failure:** Implant survival is typically 96%-99%.

Implant loss may occur at 0%-3% in routine cases and 8%-22% in complex cases.

## Peri-implantitis / Infection:

Peri-implant infection rates are generally 0%-2.7% and are associated with oral hygiene. **Marginal Bone Loss:**

Marginal bone loss is typically 0.02-1.22 mm with an average of 0.4-0.7 mm.

**Soft Tissue Recession / Esthetic Complications:** > 1 mm recession may occur in up to 11% of cases; thin biotype is a key risk factor.

**Mechanical / Prosthetic Complications:** Mechanical/prosthetic complications occur at 1%—12% and are manageable with maintenance.

## Adverse Effect

Swelling, Local inflammation, Gingival injuries, Local pain, Micro hemorrhages, Aspiration or swallowing of parts during operation, Triggering of the pharyngeal reflex during the insertion of the product, Screw fracture, Screw loosening,

Component or prosthesis failure/fracture, Hypersensitivity and/or allergic reactions, Other toxicity reactions

## Transport and Storage

The product is for single use only. It should not be used again. Store the product in a dry place at room temperature. Keep away from direct sunlight.

## Hygiene and Disinfection

DENTIMPLANT (DENTCON) is supplied Sterile and for single use only. It must not be re-cleaned or sterilized. It can damage the plumbing and design features by causing the deterioration of cleaning, disinfection and sterilization.

## Sterilization

DENTIMPLANT (DENTCON) is delivered sterile. Intact packaging protects the sterilized implant from external factors and ensures sterility until the expiry date if stored properly. When the implant is removed from the sterile packaging, asepsis rules must be followed. Damaged products in sterile packaging should not be used. Previously used or non-sterile implants should not be used under any circumstances.

## Magnetic Resonance Imaging (MRI) Safety Information

- Static magnetic fields of 1.5 and 3.0 Tesla only

- Magnetic spatial dropout field of up to 970 Gauss/cm (9.7 T/m) or less

- Whole body average specific absorption rate (SAR) of 2 W/kg (Normal Use Mode)

- Under the scanning conditions specified above, our products are expected to produce a maximum temperature increase of 3°C after 15 minutes of continuous scanning.

## More Information

In case of product complaints or for more information please contact.

Summary of Safety and Clinical Performance (SSCP) will be available in the European database. This database can be accessed at: <https://ec.europa.eu/tools/eudamed> The SSCP document related to the DENTIMPLANT abutment/implant system can be accessed through this database. Reference document: SSCP01 Dental Implant System Summary of Safety and Clinical Performance

**NOTE:** You can view the electronic version of the manuals at [\[www.dentimplant.com.tr\]](http://www.dentimplant.com.tr)

To access the SSCP file: <https://www.dentimplant.com.tr/files/>  
**Manufacturer Information**

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
**Company Address:** Ostim OSB Mah. 1269. Cad. No:17 Yenimahalle/ ANKARA

## Contact

In case of any negative impact or complaint, please contact us by including information about the complaint;  
Contact Person PRRC : Seçil KOCABIYIK

0312 424 10 28

Website: [www.dentimplant.com.tr](http://www.dentimplant.com.tr)

 Manufacturer Identification	 Production Date
 Expiry Date	 Lot number
 Reference Number	 Sterilization by Gamma Irradiation
 Do not Resterilize	 Do not Use if Package Is Damaged
 Keep Away from the Sun	 For Single Use Only
 Caution	 Check the User Manual
 MR Conditional	 CE 2292 Conformity Mark for Medical Device Regulation 2017/745
 Storage Conditions Temperature Limits	 Storage Conditions Humidity Limits
 Single sterile barrier system with protective packaging inside	 UDI Unique Device Identification
 MD Medical Device	

for e-IFU