Abutment Surgical Use Guide

Abutment

- 1- Take an impression after soft tissue is healed.
 - 2- Choose the right type and size of abutment by comparing the radiography with the taken impression and by considering adjoining teeth.
 - 3- Connect the chosen abutment to the fixture in the mouth.
 - 4- Take an impression at either fixture level or at abutment level in order to produce the final prosthesis.
 - 5- In the case of cement-retained prosthesis, use a hard finishing plaster to do a wax-up on the abutment, and manufacture the prosthesis.
 - 6- Attach the final prosthesis to the mouth and the abutment in the mouth and connect it to the fixture and then check the occlusion between the upper and lower jaws.
- 2- Uncooperative patients or patients with mental or physical disorder which may result in unstable implant, failure to fix the implant or complications in management after operation;
 - 3- Patients with the problem in bone remodeling, microcirculation or blood clotting function;
 - 4- Pregnant patients; and
 - 5- Patients with high blood pressure or diabetes mellitus.
 - 6- Patients with bone disease (osteoporosis, osteomalacia) or Bone metabolism disorders.
 - 7- Patients allergic to Titanium

3. Things to Pay Attention

1. It should be done by specialist physicians because bone tissue or transplant may be damaged due to the wrong surgical method.
2. This product should not be reused and should be applied for its original purpose.
3. Misused or damaged product should be discarded.

4. If the product is faulty, it must be sent back.

5. It should be used carefully to avoid damage and deformation of the product.

4. Side Effects

1. Delayed healing, opening, edema, anesthesia
2. Hematoma, infection, inflammation
3. May overreact.

Warnings

This product is not sterile and must be sterilized before placing in the patient's mouth. When the dentist places the final denture in the patient's mouth, he or she should check the osseointegration with x-rays, a stability checker, and percussion tests. Abutment screws should be tightened to the specified torque values. If the screw is tightened at a value within this torque range, this may damage the screw and its seat, causing it to break. Tightening at a low torque can also damage it. When connecting the screw, be especially careful not to use too much force until you are sure that the screw is seated correctly. If the initial placement of the screw is improper, it may result in damage to the screw, the abutment, or the entire implant structure. The chewing force of the prosthesis should not be applied to one side. Otherwise, stress fractures may occur.

Storage Conditions

1. Store at room temperature.
2. Keep the product away from direct sunlight and heat.
3. Store the product in a dry place.

Sterilization

Abutments are supplied to the market as non-sterile. Therefore, before abutments are used, sterilize in an autoclave at 121 °C for 15 minutes and wait until it dries.

NOTE: Using from www.dentimplant.com.tr
You can review the digital version of the manuals.

MORE INFORMATION

In case of product complaints or for more information please contact.

To access the SSCP file:
<https://www.dentimplant.com.tr/files/>

 Manufacturer Identification	 Production Date
 Expiry Date	 Lot Number
 Reference Number	 Non-Sterile
 Do not Sterilize	 Do not use if Package is Damaged
 Keep away from the Sun	 For Single Use Only
 Caution	 Check the User Manual
 Suitable for use in MR Environment.	 CE 2006 Conform to Mark for Medical Device Regulation 2017/745
 Storage Conditions Temperature Limits	 Storage Conditions Humidity Limits

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